

## Refractive Index as a PAT for Peptide and Oligonucleotide Synthesis

Beatriz G. de la Torre<sup>1</sup>, John Lopez<sup>2</sup>

<sup>1</sup> School of Laboratory Medicine and Medical Sciences, College of Health Sciences, University of KwaZulu-Natal, Durban 4041, South Africa.

<sup>2</sup> Novartis Pharma AG, Lichtstrasse 35, 4056 Basel, Switzerland

The number of peptide and oligonucleotides drugs currently on the market as well as others in advanced studies are in increasing trend. Most of them are synthesized using Solid-Phase Synthesis (SPS) technology which the main characteristic is that skips the isolation of the synthetic intermediates. The development of a real-time monitoring method would avoid the interruption of the process to take samples and run additional testing. Here we provide the demonstration that refractive index measurements of the solutions involved in SPS are a powerful Process Analytical Tool (PAT) for the real-time monitoring of the process. It will be shown how this monitoring capacity helps to the determination of the endpoint of the reactions and the optimization of all synthetic steps, which would have a direct impact on the consumption of reagents, solvents, and time, thereby contributing to greener SPS.